

产品名称: RAB23 Mouse Monoclonal Antibody
产品货号: AMM86080



产品概述 (Summary)

产品名称 (Production Name)	RAB23 Mouse Monoclonal Antibody
描述 (Description)	Mouse monoclonal Antibody
宿主 (Host)	Mouse
应用 (Application)	WB, ICC, FC
种属反应性 (Reactivity)	Human, Mouse

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	Mouse IgG1
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in TBS with 0.05% sodium azide.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	RAB23
别名 (Alternative Names)	Ras-related protein Rab-23, RAB23
基因 ID (Gene ID)	51715.0
蛋白 ID (SwissProt ID)	Q9ULC3. Purified His-tagged RAB23 protein(Fragment) was used to produced this monoclonal antibody.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:2000, ICC 1:25-1:50, FC 1:25-1:50
蛋白分子量 (Molecular Weight)	26.7kDa

研究背景 (Background)

The protein encoded by this gene belongs to the small


产品名称: RAB23 Mouse Monoclonal Antibody
产品货号: AMM86080



GTPase superfamily, Rab family. It may be involved in small GTPase mediated signal transduction and intracellular protein transportation. Alternative splicing occurs at this locus and two transcript variants encoding the same protein have been identified.

研究领域 (Research Area)

图片 (Image Data)



All lanes : Anti-RAB23 Antibody at 1:2000 dilution

注意事项 (Note)

For research use only .