

产品名称: C3 Mouse Monoclonal Antibody
产品货号: AMM86081



产品概述 (Summary)

产品名称 (Production Name)	C3 Mouse Monoclonal Antibody
描述 (Description)	Mouse monoclonal Antibody
宿主 (Host)	Mouse
应用 (Application)	WB
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	Mouse IgG1
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	C3 Complement C3, C3 and PZP-like alpha-2-macroglobulin domain-containing protein 1, Complement C3 beta chain, C3-beta-c, C3bc, Complement C3 alpha chain, C3a anaphylatoxin, Acylation stimulating protein, ASP, C3adesArg, Complement C3b alpha' chain, Complement C3c alpha' chain fragment 1, Complement C3dg fragment, Complement C3g fragment, Complement C3d fragment, Complement C3f fragment, Complement C3c alpha' chain fragment 2, C3, CPAMD1
别名 (Alternative Names)	
基因 ID (Gene ID)	718.0 P01024. This C3 antibody is generated from mice immunized with a KLH conjugated synthetic peptide between 1353-1382 amino acids from the C-terminal region of human C3.
蛋白 ID (SwissProt ID)	

产品名称: C3 Mouse Monoclonal Antibody
产品货号: AMM86081



产品应用 (Application)

稀释比 (Dilution Ratio) WB 1:1000-1:2000


蛋白分子量 (Molecular Weight) 187.1kDa

研究背景 (Background)

Complement component C3 plays a central role in the activation of complement system. Its activation is required for both classical and alternative complement activation pathways. People with C3 deficiency are susceptible to bacterial infection.

研究领域 (Research Area)

图片 (Image Data)



All lanes : Anti-C3 Antibody (C-term) at 1:2000 dilution

注意事项 (Note)

For research use only .