

产品名称: Neurogenin3 Mouse Monoclonal Antibody
产品货号: AMM86085



产品概述 (Summary)

产品名称 (Production Name)	Neurogenin3 Mouse Monoclonal Antibody
描述 (Description)	Mouse monoclonal Antibody
宿主 (Host)	Mouse
应用 (Application)	WB
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	Mouse IgG2a
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	Neurogenin3
别名 (Alternative Names)	Neurogenin-3, NGN-3, Class A basic helix-loop-helix protein 7, bHLHa7, Protein atonal homolog 5, NEUROG3, ATOH5, BHLHA7, NGN3
基因 ID (Gene ID)	50674.0
蛋白 ID (SwissProt ID)	Q9Y4Z2. This Neurogenin3 antibody is generated from a mouse immunized with a recombinant protein between 1-214 amino acids from human Neurogenin3.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:2000
蛋白分子量 (Molecular Weight)	23.1kDa

产品名称: Neurogenin3 Mouse Monoclonal Antibody
产品货号: AMM86085




研究背景 (Background)

Acts as a transcriptional regulator. Together with NKX2- 2, initiates transcriptional activation of NEUROD1. Involved in neurogenesis. Also required for the specification of a common precursor of the 4 pancreatic endocrine cell types (By similarity).

研究领域 (Research Area)

图片 (Image Data)



Anti-Neurogenin3 Antibody at 1:4000 dilution + THP-1 whole cell lysate

注意事项 (Note)

For research use only .