

产品名称: ABCB4 Mouse Monoclonal Antibody
产品货号: AMM86089



产品概述 (Summary)

产品名称 (Production Name)	ABCB4 Mouse Monoclonal Antibody
描述 (Description)	Mouse monoclonal Antibody
宿主 (Host)	Mouse
应用 (Application)	WB
种属反应性 (Reactivity)	Human, Mouse, Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	Mouse IgG1
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	ABCB4
别名 (Alternative Names)	Multidrug resistance protein 3, 3.6.3.44, ATP-binding cassette sub-family B member 4, P-glycoprotein 3, ABCB4, MDR3, PGY3
基因 ID (Gene ID)	5244.0
蛋白 ID (SwissProt ID)	P21439. This ABCB4 antibody is generated from a mouse immunized with a recombinant protein of human ABCB4.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:2000
蛋白分子量 (Molecular Weight)	141.5kDa

研究背景 (Background)


产品名称: ABCB4 Mouse Monoclonal Antibody
产品货号: AMM86089



Mediates ATP-dependent export of organic anions and drugs from the cytoplasm. Hydrolyzes ATP with low efficiency. Not capable of conferring drug resistance. Mediates the translocation of phosphatidylcholine across the canalicular membrane of the hepatocyte.

研究领域 (Research Area)

图片 (Image Data)



Anti-ABCB4 Antibody at 1:4000 dilution + MCF-7 whole cell lysate

注意事项 (Note)

For research use only .