

产品名称: DRP-2 Mouse Monoclonal Antibody
产品货号: AMM86096



产品概述 (Summary)

产品名称 (Production Name)	DRP-2 Mouse Monoclonal Antibody
描述 (Description)	Mouse monoclonal Antibody
宿主 (Host)	Mouse
应用 (Application)	WB
种属反应性 (Reactivity)	Human, Mouse, Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	Mouse IgG2b
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in TBS with 0.05% sodium azide.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	DRP-2 Dihydropyrimidinase-related protein 2, DRP-2, Collapsin response mediator protein 2, CRMP-2, N2A3, Unc-33-like phosphoprotein 2, ULIP-2, DPYSL2, CRMP2, ULIP2
别名 (Alternative Names)	
基因 ID (Gene ID)	1808.0
蛋白 ID (SwissProt ID)	Q16555. This DRP-2 antibody is generated from mice immunized with a KLH conjugated synthetic peptide between 463-492 amino acids from human DRP-2.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:2000-1:8000
蛋白分子量 (Molecular Weight)	62.3kDa


产品名称: DRP-2 Mouse Monoclonal Antibody
产品货号: AMM86096



研究背景 (Background)

研究领域 (Research Area)

图片 (Image Data)



All lanes : Anti-DRP-2 Antibody at 1:8000 dilution

注意事项 (Note)

For research use only .