

产品名称: DKK2 Mouse Monoclonal Antibody
产品货号: AMM86101



产品概述 (Summary)

产品名称 (Production Name)	DKK2 Mouse Monoclonal Antibody
描述 (Description)	Mouse monoclonal Antibody
宿主 (Host)	Mouse
应用 (Application)	WB
种属反应性 (Reactivity)	Human, Mouse

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	Mouse IgG1
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	DKK2
别名 (Alternative Names)	Dickkopf-related protein 2, Dickkopf-2, Dkk-2, hDkk-2, DKK2
基因 ID (Gene ID)	27123.0
蛋白 ID (SwissProt ID)	Q9UBU2. This DKK2 antibody is generated from a mouse immunized with a recombinant protein of human DKK2.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:2000
蛋白分子量 (Molecular Weight)	28.4kDa

研究背景 (Background)

Antagonizes canonical Wnt signaling by inhibiting LRP5/6 interaction with Wnt and by forming a ternary complex with the


产品名称: DKK2 Mouse Monoclonal Antibody
产品货号: AMM86101



transmembrane protein KREMEN that promotes internalization of LRP5/6. DKKs play an important role in vertebrate development, where they locally inhibit Wnt regulated processes such as antero- posterior axial patterning, limb development, somitogenesis and eye formation. In the adult, Dkks are implicated in bone formation and bone disease, cancer and Alzheimer disease (By similarity).

研究领域 (Research Area)

图片 (Image Data)



All lanes : Anti-DKK2 Antibody at 1:2000 dilution

注意事项 (Note)

For research use only .