

产品名称: human CD25 Mouse Monoclonal Antibody
产品货号: AMM86105



产品概述 (Summary)

产品名称 (Production Name)	human CD25 Mouse Monoclonal Antibody
描述 (Description)	Mouse monoclonal Antibody
宿主 (Host)	Mouse
应用 (Application)	FC
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	Mouse IgG2a
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in TBS with 0.05% sodium azide.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	human CD25
别名 (Alternative Names)	human CD25
基因 ID (Gene ID)	
蛋白 ID (SwissProt ID)	unknown.

产品应用 (Application)

稀释比 (Dilution Ratio)	FC 1:50-1:100
蛋白分子量 (Molecular Weight)	/


研究背景 (Background)

产品名称: human CD25 Mouse Monoclonal Antibody
产品货号: AMM86105



研究领域 (Research Area)

图片 (Image Data)



Overlay histogram showing PHA-stimulated human peripheral blood lymphocytes stained with CD25 antibody (green line). The cells were incubated in 2% bovine serum albumin to block non-specific protein-protein interactions followed by the antibody (1:50 dilution) for 60min at 37°C. The secondary antibody used was Goat Anti-Mouse IgG, DyLight® 488 Conjugated Highly Cross-Adsorbed (OJ192088) at 1/200 dilution for 40min at 37°C. Isotype control antibody (blue line) was mouse IgG2a (1 µg/1x10⁶ cells) used under the same conditions. Acquisition of >10,000 events was performed.

注意事项 (Note)

For research use only.