

产品名称: FOLR2 Mouse Monoclonal Antibody
产品货号: AMM86106



产品概述 (Summary)

产品名称 (Production Name)	FOLR2 Mouse Monoclonal Antibody
描述 (Description)	Mouse monoclonal Antibody
宿主 (Host)	Mouse
应用 (Application)	WB
种属反应性 (Reactivity)	Human, Mouse

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	Mouse IgG2b
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	FOLR2
别名 (Alternative Names)	Folate receptor beta, FR-beta, Folate receptor 2, Folate receptor, fetal/placental, Placental folate-binding protein, FBP, FOLR2
基因 ID (Gene ID)	2350.0
蛋白 ID (SwissProt ID)	P14207. This FOLR2 antibody is generated from a mouse immunized with a recombinant protein from human FOLR2.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:2000
蛋白分子量 (Molecular Weight)	29.3kDa

研究背景 (Background)


产品名称: FOLR2 Mouse Monoclonal Antibody
产品货号: AMM86106



Binds to folate and reduced folic acid derivatives and mediates delivery of 5-methyltetrahydrofolate and folate analogs into the interior of cells. Has high affinity for folate and folic acid analogs at neutral pH. Exposure to slightly acidic pH after receptor endocytosis triggers a conformation change that strongly reduces its affinity for folates and mediates their release.

研究领域 (Research Area)

图片 (Image Data)



All lanes : Anti-FOLR2 Antibody at 1:1000-1:2000 dilution

注意事项 (Note)

For research use only .