

产品名称: PRDM12 Mouse Monoclonal Antibody
产品货号: AMM86111



产品概述 (Summary)

产品名称 (Production Name)	PRDM12 Mouse Monoclonal Antibody
描述 (Description)	Mouse monoclonal Antibody
宿主 (Host)	Mouse
应用 (Application)	WB
种属反应性 (Reactivity)	Human, Mouse

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	Mouse IgG1
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in TBS with 0.05% sodium azide.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	PRDM12
别名 (Alternative Names)	PR domain zinc finger protein 12, 211-, PR domain-containing protein 12, PRDM12, PFM9
基因 ID (Gene ID)	59335.0
蛋白 ID (SwissProt ID)	Q9H4Q4. This PRDM12 antibody was raised using purified recombinant GST fusion protein encoding human PRDM12.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000
蛋白分子量 (Molecular Weight)	40.4kDa


研究背景 (Background)

产品名称: PRDM12 Mouse Monoclonal Antibody
产品货号: AMM86111



研究领域 (Research Area)

图片 (Image Data)



Anti-PRDM12 Antibody at 1:1000 dilution + 293 whole cell lysate

注意事项 (Note)

For research use only .