

产品名称: ENOG Mouse Monoclonal Antibody
产品货号: AMM86114



产品概述 (Summary)

产品名称 (Production Name)	ENOG Mouse Monoclonal Antibody
描述 (Description)	Mouse monoclonal Antibody
宿主 (Host)	Mouse
应用 (Application)	WB,IHC
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	Mouse IgG2b
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in TBS with 0.05% sodium azide.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	ENOG
别名 (Alternative Names)	Gamma-enolase, 2-phospho-D-glycerate hydro-lyase, Enolase 2, Neural enolase, Neuron-specific enolase, NSE, ENO2
基因 ID (Gene ID)	2026.0
蛋白 ID (SwissProt ID)	P09104.This ENOG monoclonal antibody is generated from mouse immunized with ENOG recombinant protein.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000,IHC 1:100-1:200
蛋白分子量 (Molecular Weight)	47.3kDa

研究背景 (Background)


产品名称: ENOG Mouse Monoclonal Antibody
产品货号: AMM86114



This gene encodes one of the three enolase isoenzymes found in mammals. This isoenzyme, a homodimer, is found in mature neurons and cells of neuronal origin. A switch from alpha enolase to gamma enolase occurs in neural tissue during development in rats and primates.

研究领域 (Research Area)

图片 (Image Data)



All lanes : Anti-ENOG at 1:1000 dilution

注意事项 (Note)

For research use only .