

产品名称: GDE1 Mouse Monoclonal Antibody
产品货号: AMM86116



产品概述 (Summary)

产品名称 (Production Name)	GDE1 Mouse Monoclonal Antibody
描述 (Description)	Mouse monoclonal Antibody
宿主 (Host)	Mouse
应用 (Application)	WB
种属反应性 (Reactivity)	Human, Mouse

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	Mouse IgG1
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	GDE1
别名 (Alternative Names)	Glycerophosphodiester phosphodiesterase 1, 3.1.4.44, Membrane-interacting protein of RGS16, RGS16-interacting membrane protein, GDE1, MIR16
基因 ID (Gene ID)	51573.0
蛋白 ID (SwissProt ID)	Q9NZC3. This GDE1 antibody is generated from a mouse immunized with a recombinant protein from the human region of human GDE1.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:2000
蛋白分子量 (Molecular Weight)	37.7kDa

研究背景 (Background)


产品名称: GDE1 Mouse Monoclonal Antibody
产品货号: AMM86116



Has glycerophosphoinositol phosphodiesterase activity. Has little or no activity towards glycerophosphocholine. GDE1 activity can be modulated by G-protein signaling pathways (By similarity).

研究领域 (Research Area)

图片 (Image Data)



Anti-GDE1 Antibody at 1:2000 dilution + human brain tissue lysate

注意事项 (Note)

For research use only .