

产品名称: MAP3K13 (LZK) Mouse Monoclonal Antibody
产品货号: AMM86117



产品概述 (Summary)

产品名称 (Production Name)	MAP3K13 (LZK) Mouse Monoclonal Antibody
描述 (Description)	Mouse monoclonal Antibody
宿主 (Host)	Mouse
应用 (Application)	WB
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	Mouse IgM
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	MAP3K13 (LZK)
别名 (Alternative Names)	Mitogen-activated protein kinase kinase kinase 13, Leucine zipper-bearing kinase, Mixed lineage kinase, MLK, MAP3K13 (http://www.genenames.org/cgi-bin/gene_symbol_report?hgnc_id=6852 target="_blank">HGNC:6852)
基因 ID (Gene ID)	9175.0
蛋白 ID (SwissProt ID)	O43283. Purified His-tagged MAP3K13 protein fragment was used to produced this monoclonal antibody.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000
蛋白分子量 (Molecular Weight)	108.3kDa

产品名称: MAP3K13 (LZK) Mouse Monoclonal Antibody
产品货号: AMM86117




研究背景 (Background)

Activates the JUN N-terminal pathway through activation of the MAP kinase kinase MAP2K7. Acts synergistically with PRDX3 to regulate the activation of NF-kappa-B in the cytosol. This activation is kinase-dependent and involves activating the IKK complex, the IKBKB-containing complex that phosphorylates inhibitors of NF-kappa-B.

研究领域 (Research Area)

MAPK signaling pathway

图片 (Image Data)



Anti-MAP3K13 (LZK) Antibody (C-term) at 1:500 dilution + NCI-H460 whole cell lysate

注意事项 (Note)

For research use only .