

产品名称: ITGA11 Mouse Monoclonal Antibody
产品货号: AMM86119



产品概述 (Summary)

产品名称 (Production Name)	ITGA11 Mouse Monoclonal Antibody
描述 (Description)	Mouse monoclonal Antibody
宿主 (Host)	Mouse
应用 (Application)	WB
种属反应性 (Reactivity)	Human, Mouse

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	Mouse IgG2b
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	ITGA11
别名 (Alternative Names)	Integrin alpha-11, ITGA11
基因 ID (Gene ID)	22801.0
蛋白 ID (SwissProt ID)	Q9UKX5. This ITGA11 antibody is generated from a mouse immunized with a KLH conjugated synthetic peptide between 675-880 amino acids from human ITGA11.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:2000-1:8000
蛋白分子量 (Molecular Weight)	133.5kDa

研究背景 (Background)

产品名称: ITGA11 Mouse Monoclonal Antibody
产品货号: AMM86119




Integrin alpha-11/beta-1 is a receptor for collagen.

研究领域 (Research Area)

PI3K-Akt signaling pathway

图片 (Image Data)



All lanes : Anti-ITGA11 Antibody at 1:8000 dilution

注意事项 (Note)

For research use only .