

产品名称: GSTA1 Mouse Monoclonal Antibody
产品货号: AMM86122



产品概述 (Summary)

产品名称 (Production Name)	GSTA1 Mouse Monoclonal Antibody
描述 (Description)	Mouse monoclonal Antibody
宿主 (Host)	Mouse
应用 (Application)	WB
种属反应性 (Reactivity)	Human, Mouse, Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	Mouse IgG1
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	GSTA1
别名 (Alternative Names)	Glutathione S-transferase A1, 2.5.1.18, GST HA subunit 1, GST class-alpha member 1, GST-epsilon, GSTA1-1, GTH1, Glutathione S-transferase A1, N-terminally processed, GSTA1
基因 ID (Gene ID)	2938.0
蛋白 ID (SwissProt ID)	P08263. This GSTA1 antibody is generated from a mouse immunized with a recombinant protein from the human region of human GSTA1.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:2000-1:5000
蛋白分子量 (Molecular Weight)	25.6kDa

产品名称: GSTA1 Mouse Monoclonal Antibody
产品货号: AMM86122




研究背景 (Background)

Conjugation of reduced glutathione to a wide number of exogenous and endogenous hydrophobic electrophiles.

研究领域 (Research Area)

图片 (Image Data)



All lanes : Anti-GSTA1 Antibody at 1:5000 dilution

注意事项 (Note)

For research use only .