

产品名称: VWF Mouse Monoclonal Antibody  
产品货号: AMM86125



## 产品概述 (Summary)

|                        |                               |
|------------------------|-------------------------------|
| 产品名称 (Production Name) | VWF Mouse Monoclonal Antibody |
| 描述 (Description)       | Mouse monoclonal Antibody     |
| 宿主 (Host)              | Mouse                         |
| 应用 (Application)       | WB                            |
| 种属反应性 (Reactivity)     | Human                         |

## 产品性能 (Performance)

|                     |                                                                                          |
|---------------------|------------------------------------------------------------------------------------------|
| 偶联物 (Conjugation)   | Unconjugated                                                                             |
| 修饰 (Modification)   | Unmodified                                                                               |
| 同种型 (Isotype)       | Mouse IgG1                                                                               |
| 克隆 (Clonality)      | Monoclonal                                                                               |
| 形式 (Form)           | Liquid                                                                                   |
| 存放说明 (Storage)      | Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles. |
| 储存溶液 (Buffer)       | Purified antibody in PBS with 0.05% sodium azide.                                        |
| 纯化方式 (Purification) | Affinity Purification                                                                    |

## 免疫原信息 (Immunogen)

|                        |                                                                                             |
|------------------------|---------------------------------------------------------------------------------------------|
| 基因名 (Gene Name)        | VWF                                                                                         |
| 别名 (Alternative Names) | von Willebrand factor, vWF, von Willebrand antigen 2, von Willebrand antigen II, VWF, F8VWF |
| 基因 ID (Gene ID)        | 7450.0                                                                                      |
| 蛋白 ID (SwissProt ID)   | P04275. This antibody is generated from a mouse immunized with .VWF recombinant protein.    |

## 产品应用 (Application)

|                          |                  |
|--------------------------|------------------|
| 稀释比 (Dilution Ratio)     | WB 1:1000-1:2000 |
| 蛋白分子量 (Molecular Weight) | 309.3kDa         |

## 研究背景 (Background)


**产品名称: VWF Mouse Monoclonal Antibody**  
**产品货号: AMM86125**



Important in the maintenance of hemostasis, it promotes adhesion of platelets to the sites of vascular injury by forming a molecular bridge between sub-endothelial collagen matrix and platelet-surface receptor complex GPIb-IX-V. Also acts as a chaperone for coagulation factor VIII, delivering it to the site of injury, stabilizing its heterodimeric structure and protecting it from premature clearance from plasma.

## 研究领域 (Research Area)

## 图片 (Image Data)



Anti-VWF at dilution + VWF whole cell lysate

## 注意事项 (Note)

For research use only .