

产品名称: SET07 Mouse Monoclonal Antibody
产品货号: AMM86128



产品概述 (Summary)

产品名称 (Production Name)	SET07 Mouse Monoclonal Antibody
描述 (Description)	Mouse monoclonal Antibody
宿主 (Host)	Mouse
应用 (Application)	WB,IHC,ICC
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	Mouse IgG1
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	SET07 N-lysine methyltransferase SETD8, 211-, H4-K20-HMTase SETD8, Histone-lysine N-methyltransferase SETD8, Lysine N-methyltransferase 5A, PR/SET domain-containing protein 07, PR-Set7, PR/SET07, SET domain-containing protein 8, SETD8, KMT5A, PRSET7, SET07, SET8
别名 (Alternative Names)	
基因 ID (Gene ID)	387893.0
蛋白 ID (SwissProt ID)	Q9NQR1.This SET07 antibody was raised using purified recombinant GST fusion protein encoding the N-terminal region of human SET07.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:2000,IHC 1:100-1:500,ICC 1:20-1:50
蛋白分子量 (Molecular Weight)	42.9kDa


产品名称: SET07 Mouse Monoclonal Antibody
产品货号: AMM86128



研究背景 (Background)

研究领域 (Research Area)

图片 (Image Data)



All lanes :SET07 Antibody at 1:2000 dilution

注意事项 (Note)

For research use only .