

产品名称: PSMC2 Mouse Monoclonal Antibody
产品货号: AMM86136



产品概述 (Summary)

产品名称 (Production Name)	PSMC2 Mouse Monoclonal Antibody
描述 (Description)	Mouse monoclonal Antibody
宿主 (Host)	Mouse
应用 (Application)	WB
种属反应性 (Reactivity)	Human, Mouse

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	Mouse IgG1
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	PSMC2
别名 (Alternative Names)	26S protease regulatory subunit 7, 26S proteasome AAA-ATPase subunit RPT1, Proteasome 26S subunit ATPase 2, Protein MSS1, PSMC2, MSS1
基因 ID (Gene ID)	5701.0
蛋白 ID (SwissProt ID)	P35998. This antibody is generated from a mouse immunized with a recombinant protein from human.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:2000
蛋白分子量 (Molecular Weight)	48.6kDa

研究背景 (Background)


产品名称: PSMC2 Mouse Monoclonal Antibody
产品货号: AMM86136



The 26S protease is involved in the ATP-dependent degradation of ubiquitinated proteins. The regulatory (or ATPase) complex confers ATP dependency and substrate specificity to the 26S complex. In case of HIV-1 infection, positive modulator of Tat-mediated transactivation.

研究领域 (Research Area)

图片 (Image Data)



All lanes : Anti-PSMC2 Antibody at 1:1000 dilution + C6 cell lysate

注意事项 (Note)

For research use only .