

产品名称: MTHFD2 Mouse Monoclonal Antibody
产品货号: AMM86137



产品概述 (Summary)

产品名称 (Production Name)	MTHFD2 Mouse Monoclonal Antibody
描述 (Description)	Mouse monoclonal Antibody
宿主 (Host)	Mouse
应用 (Application)	WB
种属反应性 (Reactivity)	Human, Mouse

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	Mouse IgG1
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	MTHFD2
别名 (Alternative Names)	Bifunctional methylenetetrahydrofolate dehydrogenase/cyclohydrolase, mitochondrial, NAD-dependent methylenetetrahydrofolate dehydrogenase, 1.5.1.15, Methenyltetrahydrofolate cyclohydrolase, 3.5.4.9, MTHFD2, NMDMC
基因 ID (Gene ID)	10797.0
蛋白 ID (SwissProt ID)	P13995. This MTHFD2 antibody is generated from a mouse immunized with a KLH conjugated synthetic peptide between 266-299 amino acids from the C-terminal region of human MTHFD2.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:2000
蛋白分子量 (Molecular Weight)	37.9kDa


产品名称: MTHFD2 Mouse Monoclonal Antibody
产品货号: AMM86137



研究背景 (Background)

研究领域 (Research Area)

图片 (Image Data)



All lanes : Anti-MTHFD2 Antibody (C-term) at 1:2000 dilution

注意事项 (Note)

For research use only .