

产品名称: Phosphoserine(5B12) Mouse Monoclonal Antibody
产品货号: AMM86138



产品概述 (Summary)

产品名称 (Production Name)	Phosphoserine(5B12) Mouse Monoclonal Antibody
描述 (Description)	Mouse monoclonal Antibody
宿主 (Host)	Mouse
应用 (Application)	WB,IHC
种属反应性 (Reactivity)	Other

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Phosphorylated
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide, 0.5% protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	Phosphoserine(5B12)
别名 (Alternative Names)	Phosphoserine(5B12)
基因 ID (Gene ID)	
蛋白 ID (SwissProt ID)	/Recombinant Protein

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:2000-1:5000, IHC 1:50-1:200
蛋白分子量 (Molecular Weight)	/


研究背景 (Background)

产品名称: Phosphoserine(5B12) Mouse Monoclonal Antibody
产品货号: AMM86138




研究领域 (Research Area)


图片 (Image Data)



Western blot analysis of 1) Hela, 2) Rat Brain Tissue, 3) Mouse Brain Tissue with Phosphoserine Mouse mAb diluted at 1:2,000.



Immunohistochemical analysis of paraffin-embedded Human Lung Carcinoma using Phosphoserine Mouse mAb diluted at 1:200.



Immunohistochemical analysis of paraffin-embedded Human Breast Carcinoma using Phosphoserine Mouse mAb diluted at 1:200.

产品名称: Phosphoserine(5B12) Mouse Monoclonal
Antibody
产品货号: AMM86138



1:200.

注意事项 (Note)

For research use only .