

产品名称: Phospho-Akt (S473)Mouse Monoclonal Antibody
产品货号: AMM86142



产品概述 (Summary)

产品名称 (Production Name)	Phospho-Akt (S473)Mouse Monoclonal Antibody
描述 (Description)	Mouse monoclonal Antibody
宿主 (Host)	Mouse
应用 (Application)	WB,IHC
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Phosphorylated
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide,0.5%protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	Phospho-Akt (S473)
别名 (Alternative Names)	AKT1
基因 ID (Gene ID)	207;
蛋白 ID (SwissProt ID)	P31749.Synthetic Peptide

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:2000,IHC 1:100-1:200
蛋白分子量 (Molecular Weight)	60kDa

研究背景 (Background)


产品名称: Phospho-Akt (S473) Mouse Monoclonal Antibody
产品货号: AMM86142




研究领域 (Research Area)

Apoptosis, TGF-beta signaling pathway, PI3K-Akt signaling pathway, mTOR signaling pathway, MAPK signaling pathway, Jak-STAT signaling pathway

图片 (Image Data)



Western blot analysis of PC3 Cell Lysate using Phospho-Akt Ser473 Mouse mAb diluted at 1:1000.



Immunohistochemical analysis of paraffin-embedded Human Lung Carcinoma Tissue using Phospho-Akt Ser473 Mouse mAb diluted at 1:200.

注意事项 (Note)

For research use only .