

产品名称: CK17 Mouse Monoclonal Antibody
产品货号: AMM86154



产品概述 (Summary)

产品名称 (Production Name)	CK17 Mouse Monoclonal Antibody
描述 (Description)	Mouse monoclonal Antibody
宿主 (Host)	Mouse
应用 (Application)	WB,IHC,ICC,IP
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	Mouse IgG1
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	CK17
别名 (Alternative Names)	CK17 Mouse
基因 ID (Gene ID)	
蛋白 ID (SwissProt ID)	Q04695.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:2000,IHC 1:50-1:300,ICC 1:100-1:200,IP 1:20-1:300
蛋白分子量 (Molecular Weight)	48kDa

研究背景 (Background)

Cytokeratin17 (CK17) is a protein that in humans is encoded by the KRT17 gene. It is found in nail beds, hair follicles, sebaceous glands, and other epidermal appendages. Mutations in the gene encoding this protein lead to Jackson-Lawler


产品名称: CK17 Mouse Monoclonal Antibody
产品货号: AMM86154




type pachyonychia congenita and steatocystoma multiplex.

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of 1) HeLa, 2) MCF7, 3) 293T, using AMM86154 diluted at 1:2,000.



1 Input: HeLa Cell Lysate 2 IP product: IP dilute 1:200 Western blot analysis: primary antibody : AMM86154 1:1,000 Secondary antibody: Goat anti-Mouse IgG, Light chain specific, 1:5,000

注意事项 (Note)

For research use only .