

产品名称: Histone H2B Mouse Monoclonal Antibody
产品货号: AMM86158



产品概述 (Summary)

产品名称 (Production Name)	Histone H2B Mouse Monoclonal Antibody
描述 (Description)	Mouse monoclonal Antibody
宿主 (Host)	Mouse
应用 (Application)	WB, ICC
种属反应性 (Reactivity)	Human, Mouse, Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	Histone H2B
别名 (Alternative Names)	H2B GL105 antibody, H2B.1 antibody, Histone H2B type 1B antibody, HIST1H2BB antibody
基因 ID (Gene ID)	3018;
蛋白 ID (SwissProt ID)	P33778.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:3000, ICC 1:100-1:200
蛋白分子量 (Molecular Weight)	14kDa

研究背景 (Background)


Histone H2B is one of the 5 main histone proteins involved in the structure of chromatin in eukaryotic cells.

产品名称: Histone H2B Mouse Monoclonal Antibody
产品货号: AMM86158



研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of 1) Hela, 2) 3T3, 3) Raw264.7, 4) Rat Brain, 5) Rat Kidney, using AMM86158 diluted at 1:2,000.

注意事项 (Note)

For research use only .