

产品名称: HAO1 Mouse Monoclonal Antibody
产品货号: AMM86161



产品概述 (Summary)

产品名称 (Production Name)	HAO1 Mouse Monoclonal Antibody
描述 (Description)	Mouse monoclonal Antibody
宿主 (Host)	Mouse
应用 (Application)	WB
种属反应性 (Reactivity)	Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide,0.5%protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	HAO1
别名 (Alternative Names)	HAO1
基因 ID (Gene ID)	54363;
蛋白 ID (SwissProt ID)	Q9UJM8.Recombinant Protein

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000
蛋白分子量 (Molecular Weight)	41kDa


研究背景 (Background)

产品名称: HAO1 Mouse Monoclonal Antibody
产品货号: AMM86161



研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of 1) Mouse Liver Tissue, 2) Rat Liver Tissue using HAO1 Monoclonal Antibody.

注意事项 (Note)

For research use only .