

产品名称: Pan Methylated Lysine Mouse Monoclonal Antibody
产品货号: AMM86167



产品概述 (Summary)

产品名称 (Production Name)	Pan Methylated Lysine Mouse Monoclonal Antibody
描述 (Description)	Mouse monoclonal Antibody
宿主 (Host)	Mouse
应用 (Application)	WB,IHC
种属反应性 (Reactivity)	Other

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	Pan Methylated Lysine
别名 (Alternative Names)	Pan Methylated Lysine
基因 ID (Gene ID)	
蛋白 ID (SwissProt ID)	P68431.Conjugated Protein

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:2000,IHC 1:200-1:500
蛋白分子量 (Molecular Weight)	Calculated MW: 15 kDa; Observed MW: 15 kDa


研究背景 (Background)

产品名称: Pan Methylated Lysine Mouse Monoclonal Antibody
产品货号: AMM86167




研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of HeLa using Pan Methylated Lysine Monoclonal Antibody.



Immunohistochemical analysis of paraffin-embedded Human Breast Carcinoma using Pan Methylated Lysine Monoclonal Antibody.

注意事项 (Note)

For research use only .