

产品名称: Histone H3 (Tri-Methyl-Lys79) Mouse Monoclonal Antibody
产品货号: AMM86168



产品概述 (Summary)

产品名称 (Production Name)	Histone H3 (Tri-Methyl-Lys79) Mouse Monoclonal Antibody
描述 (Description)	Mouse monoclonal Antibody
宿主 (Host)	Mouse
应用 (Application)	WB,IP
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	Histone H3 (Tri-Methyl-Lys79)
别名 (Alternative Names)	HIST1H3A; H3FA; HIST1H3B; H3FL; HIST1H3C; H3FC; HIST1H3D; H3FB; HIST1H3E; H3FD; HIST1H3F; H3FI; HIST1H3G; H3FH; HIST1H3H; H3FK; HIST1H3I; H3FF; HIST1H3J; H3FJ; Histone H3.1; Histone H3/a; Histone H3/b; Histone H3/c; Histone H3/d; Histone H3/f; Histone H3/
基因 ID (Gene ID)	8350;8351;8352;8353;8354;8355;8356;8357;8358;8968;
蛋白 ID (SwissProt ID)	P68431.tri methyl peptide at K79 of Histone H3

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,IP 1:20-1:300
蛋白分子量 (Molecular Weight)	15kDa


产品名称: Histone H3 (Tri-Methyl-Lys79) Mouse Monoclonal Antibody
产品货号: AMM86168




研究背景 (Background)

研究领域 (Research Area)

图片 (Image Data)



1) Input: HeLa Cell Lysate



1) Input: HeLa Cell Lysate

注意事项 (Note)

For research use only .