

产品名称: Acetyl-NF kB P65(K314/K315) (2A11) Mouse
Monoclonal Antibody
产品货号: AMM86170



产品概述 (Summary)

产品名称 (Production Name)	Acetyl-NF kB P65(K314/K315) (2A11) Mouse Monoclonal Antibody
描述 (Description)	Mouse monoclonal Antibody
宿主 (Host)	Mouse
应用 (Application)	IHC
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Acetylated
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide,0.5%protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	Acetyl NF kB P65(K314/K315) (2A11)
别名 (Alternative Names)	Transcription factor p65 (Nuclear factor NF-kappa-B p65 subunit) (Nuclear factor of kappa light polypeptide gene enhancer in B-cells 3)
基因 ID (Gene ID)	5970;
蛋白 ID (SwissProt ID)	Q04206.Synthetic Peptide of Acetyl NF kB P65(K314/K315)

产品应用 (Application)

稀释比 (Dilution Ratio)	IHC 1:50-1:200
蛋白分子量 (Molecular Weight)	65kDa

产品名称: Acetyl-NF kB P65(K314/K315) (2A11) Mouse
Monoclonal Antibody
产品货号: AMM86170




研究背景 (Background)


研究领域 (Research Area)

PI3K-Akt signaling pathway, MAPK signaling pathway

图片 (Image Data)



Immunohistochemical analysis of paraffin-embedded Human Breast Carcinoma Tissue using Acetyl NF kB P65(K314/K315)
Mouse mAb diluted at 1:200.



Immunohistochemical analysis of paraffin-embedded Human Lung Carcinoma Tissue using Acetyl NF kB P65(K314/K315)
Mouse mAb diluted at 1:200.

注意事项 (Note)

For research use only .