

产品名称: Cbfa2t2 Rat Monoclonal Antibody
产品货号: AMR83163



产品概述 (Summary)

产品名称 (Production Name)	Cbfa2t2 Rat Monoclonal Antibody
描述 (Description)	Rat monoclonal antibody
宿主 (Host)	Rat
应用 (Application)	IHC
种属反应性 (Reactivity)	Mouse, Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	Rat IgG2c
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	Cbfa2t2
别名 (Alternative Names)	Cbfa2t2
基因 ID (Gene ID)	
蛋白 ID (SwissProt ID)	O70374.MTGR1 peptide

产品应用 (Application)

稀释比 (Dilution Ratio)	IHC 1:50-1:100
蛋白分子量 (Molecular Weight)	/


研究背景 (Background)

产品名称: Cbfa2t2 Rat Monoclonal Antibody
产品货号: AMR83163



研究领域 (Research Area)

图片 (Image Data)



IHC-P analysis of testis tissue by anti-Cbfa2t2 antibody (AMR83163)

注意事项 (Note)

For research use only .