

产品名称: Ctf Rat Monoclonal Antibody
产品货号: AMR83247



产品概述 (Summary)

产品名称 (Production Name)	Ctf Rat Monoclonal Antibody
描述 (Description)	Rat monoclonal antibody
宿主 (Host)	Rat
应用 (Application)	WB,IP
种属反应性 (Reactivity)	Mouse

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	Rat IgG2b
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	Ctf
别名 (Alternative Names)	Ctf
基因 ID (Gene ID)	
蛋白 ID (SwissProt ID)	Q61164.Purified recombinant fragment of mouse Ctf

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000,IP 1:100-1:200
蛋白分子量 (Molecular Weight)	135kDa


研究背景 (Background)

产品名称: Ctf Rat Monoclonal Antibody
产品货号: AMR83247




研究领域 (Research Area)


图片 (Image Data)



Lane 1: 15 µg of HeLa lysate; Lane 2: 15 µg of Jurkat lysate; Lane 3: 15 µg of MOLT-4 lysate; Result: AMR83247 can detect CTcf by Western blotting;



Western blot analysis using CTcf mouse mAb against NIH/3T3 lysate



Lane 1: MOLT-4 lysate; Lane 2: CTcf immunoprecipitated from MOLT-4 lysate by AMR83247; Lane 3: The same as Lane 2 but use IgG isotype control antibody; Result: AMR83247 can immunoprecipitate CTcf;

产品名称: Ctf Rat Monoclonal Antibody
产品货号: AMR83247



注意事项 (Note)

For research use only .