

产品名称: Smc3 Rat Monoclonal Antibody
产品货号: AMR83258



产品概述 (Summary)

产品名称 (Production Name)	Smc3 Rat Monoclonal Antibody
描述 (Description)	Rat monoclonal antibody
宿主 (Host)	Rat
应用 (Application)	WB,IHC
种属反应性 (Reactivity)	Mouse

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	Rat IgG2a
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	Smc3
别名 (Alternative Names)	Smc3
基因 ID (Gene ID)	
蛋白 ID (SwissProt ID)	Q9CW03.Purified recombinant fragment of mouse Smc3

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000,IHC 1:100-1:200
蛋白分子量 (Molecular Weight)	55kDa


研究背景 (Background)

产品名称: Smc3 Rat Monoclonal Antibody
产品货号: AMR83258




研究领域 (Research Area)

图片 (Image Data)



Western blot analysis using SMC3 mAb against mouse testicles lysate and NIH/3T3 lysate



IHC-P analysis of colon tissue by anti-SMC3 antibody (AMR83258)

注意事项 (Note)

For research use only .