

产品名称: FGA Rat Monoclonal Antibody
产品货号: AMR83289



产品概述 (Summary)

产品名称 (Production Name)	FGA Rat Monoclonal Antibody
描述 (Description)	Rat monoclonal antibody
宿主 (Host)	Rat
应用 (Application)	WB,IHC,IP
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	Rat IgG2a
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	FGA
别名 (Alternative Names)	FGA
基因 ID (Gene ID)	
蛋白 ID (SwissProt ID)	P02671.Purified recombinant fragment of human FGA

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000,IHC 1:200-1:500,IP 1:100-1:200
蛋白分子量 (Molecular Weight)	95kDa


研究背景 (Background)

产品名称: FGA Rat Monoclonal Antibody
产品货号: AMR83289




研究领域 (Research Area)


图片 (Image Data)



Western blot analysis using human FGA mAb against human plasma



Lane 1: plasma; Lane 2: FGA immunoprecipitated from plasma by AMR83289; Lane 3: The same as Lane 2 but use IgG isotype control antibody; Result: AMR83289 can immunoprecipitate FGA;



IHC-P analysis of placenta tissue by FGA antibody (AMR83289)

注意事项 (Note)

产品名称: FGA Rat Monoclonal Antibody
产品货号: AMR83289



For research use only .