

产品名称: Mapre2 Rat Monoclonal Antibody
产品货号: AMR83412



产品概述 (Summary)

产品名称 (Production Name)	Mapre2 Rat Monoclonal Antibody
描述 (Description)	Rat monoclonal antibody
宿主 (Host)	Rat
应用 (Application)	WB,IHC,IP
种属反应性 (Reactivity)	Mouse, Human, Rat, Ape, Hamster

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	Rat IgG2b
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	Mapre2
别名 (Alternative Names)	Mapre2
基因 ID (Gene ID)	
蛋白 ID (SwissProt ID)	Q8R001.Purified recombinant fragment of mouse Mapre2

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000,IHC 1:100-1:200,IP 1:100-1:200
蛋白分子量 (Molecular Weight)	41kDa,43kDa


研究背景 (Background)

产品名称: Mapre2 Rat Monoclonal Antibody
产品货号: AMR83412




研究领域 (Research Area)


图片 (Image Data)



Lane 1: 15 µg of rat brain tissue lysate; Lane 2: 15 µg of mouse brain tissue lysate; Lane 3: 15 µg of NIH-3T3 lysate; Lane 4: 15 µg of PMA treated Jurkat lysate; Result: AMR83412 can detect Mapre2 by Western blotting.



Lane 1: Rat brain tissue lysate; Lane 2: Mapre2 immunoprecipitated from rat brain tissue lysate by AMR83412; Lane 3: The same as Lane 2 but use IgG isotype control antibody; Result: AMR83412 can immunoprecipitate Mapre2;



IHC-P analysis of mouse brain tissue by Mapre2 Antibody (AMR83412)

产品名称: Mapre2 Rat Monoclonal Antibody
产品货号: AMR83412



注意事项 (Note)

For research use only .