

产品名称: Histone H3K9ac Rat Monoclonal Antibody
产品货号: AMR83415



产品概述 (Summary)

产品名称 (Production Name)	Histone H3K9ac Rat Monoclonal Antibody
描述 (Description)	Rat monoclonal antibody
宿主 (Host)	Rat
应用 (Application)	WB,IHC,IP
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	Rat IgG1
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	Histone H3K9ac
别名 (Alternative Names)	Histone H3K9ac
基因 ID (Gene ID)	
蛋白 ID (SwissProt ID)	P68431/P84243.Purified recombinant fragment of human Histone H3K9ac

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000,IHC 1:200-1:500,IP 1:100-1:200
蛋白分子量 (Molecular Weight)	17kDa


研究背景 (Background)

产品名称: Histone H3K9ac Rat Monoclonal Antibody
产品货号: AMR83415



研究领域 (Research Area)

图片 (Image Data)



Microtiter wells were coated with various peptides. AMR83415 was used as the primary antibody and peroxidase conjugated goat anti-rat IgG was used as the the secondary antibody.

注意事项 (Note)

For research use only .