

产品名称: MAPRE3 Rat Monoclonal Antibody
产品货号: AMR83580



产品概述 (Summary)

产品名称 (Production Name)	MAPRE3 Rat Monoclonal Antibody
描述 (Description)	Rat monoclonal antibody
宿主 (Host)	Rat
应用 (Application)	IHC,IP
种属反应性 (Reactivity)	Human, Mouse

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	Rat IgG2a
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	MAPRE3
别名 (Alternative Names)	MAPRE3
基因 ID (Gene ID)	
蛋白 ID (SwissProt ID)	Q9UPY8.Purified recombinant fragment of human MAPRE3

产品应用 (Application)

稀释比 (Dilution Ratio)	IHC 1:200-1:500,IP 1:100-1:200
蛋白分子量 (Molecular Weight)	/


研究背景 (Background)

产品名称: MAPRE3 Rat Monoclonal Antibody
产品货号: AMR83580




研究领域 (Research Area)

图片 (Image Data)



Lane 1: Mouse brain lysate ; Lane 2: MAPRE3 immunoprecipitated from mouse brain lysate by AMR83580 ; Lane 3: The same Lane 2 but use IgG isotype control antibody ;Result: AMR83580 can immunoprecipitate MAPRE3;



IHC-P analysis of lung tissue by MAPRE3 antibody (AMR83580)

注意事项 (Note)

For research use only .