

产品名称: Clasp1 Rat Monoclonal Antibody
产品货号: AMR83591



产品概述 (Summary)

产品名称 (Production Name)	Clasp1 Rat Monoclonal Antibody
描述 (Description)	Rat monoclonal antibody
宿主 (Host)	Rat
应用 (Application)	WB,IP
种属反应性 (Reactivity)	Mouse, Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	Rat IgG2a
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	Clasp1
别名 (Alternative Names)	Clasp1
基因 ID (Gene ID)	
蛋白 ID (SwissProt ID)	Q80TV8.Purified recombinant fragment of mouse Clasp1

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000,IP 1:100-1:200
蛋白分子量 (Molecular Weight)	180kDa


研究背景 (Background)

产品名称: Clasp1 Rat Monoclonal Antibody
产品货号: AMR83591




研究领域 (Research Area)

图片 (Image Data)



Lane 1: 15 µg of U-937 lysate; Lane 2: 15 µg of K-562 lysate; Lane 3: 15 µg of rat heart tissue lysate; Result: AMR83591 can detect Clasp1 by Western blotting;



Lane 1: K-562 lysate ; Lane 2: Clasp1 immunoprecipitated from K-562 lysate by AMR83591 ; Lane 3: The same as Lane 2 but use IgG isotype control antibody ; Result: AMR83591 can immunoprecipitate Clasp1;

注意事项 (Note)

For research use only .