

产品名称: Tmem176b Rat Monoclonal Antibody
产品货号: AMR83673



产品概述 (Summary)

产品名称 (Production Name)	Tmem176b Rat Monoclonal Antibody
描述 (Description)	Rat monoclonal antibody
宿主 (Host)	Rat
应用 (Application)	WB
种属反应性 (Reactivity)	Mouse

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	Rat IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	Tmem176b
别名 (Alternative Names)	Tmem176b
基因 ID (Gene ID)	
蛋白 ID (SwissProt ID)	Q9R1Q6.Purified recombinant fragment of mouse Tmem176b

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:2000
蛋白分子量 (Molecular Weight)	26kDa, 28kDa


研究背景 (Background)

产品名称: Tmem176b Rat Monoclonal Antibody
产品货号: AMR83673



研究领域 (Research Area)

图片 (Image Data)



Western blot analysis using Tmem176b mAb against mouse liver tissue lysate and mouse intestine tissue lysate

注意事项 (Note)

For research use only .