

产品名称: KDEL Rabbit Polyclonal Antibody
产品货号: APRab00011



产品概述 (Summary)

产品名称 (Production Name)	KDEL Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ICC/IF,FC
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
纯化方式 (Purification)	Affinity Chromatography

免疫原信息 (Immunogen)

基因名 (Gene Name)	KDEL1
别名 (Alternative Names)	Lys-Asp-Glu-Leu (KDEL
基因 ID (Gene ID)	10945
蛋白 ID (SwissProt ID)	P24390.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000,IHC 1:50-1:100,ICC/IF 1:50-1:200,FC 1:50-1:100
蛋白分子量 (Molecular Weight)	Calculated MW: 25 kDa; Observed MW: 78kDa kDa

研究背景 (Background)

Determines the specificity of the luminal ER protein retention system.Required for the retention of luminal endoplasmic

产品名称: KDEL Rabbit Polyclonal Antibody
产品货号: APRab00011




reticulum proteins. Determines the specificity of the luminal ER protein retention system. Also required for normal vesicular traffic through the Golgi. This receptor recognizes K-D-E-L.


研究领域 (Research Area)

Tags & Cell Markers

图片 (Image Data)



Western blot analysis of KDEL in HeLa lysates using KDEL antibody.



Immunohistochemistry analysis of paraffin-embedded Human liver using KDEL antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.

注意事项 (Note)

For research use only .