

产品名称: OTX1 Rabbit Polyclonal Antibody
产品货号: APRab00018



产品概述 (Summary)

产品名称 (Production Name)	OTX1 Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
纯化方式 (Purification)	Affinity Chromatography

免疫原信息 (Immunogen)

基因名 (Gene Name)	OTX1
别名 (Alternative Names)	Homeobox protein OTX1 (Orthodenticle homolog 1)
基因 ID (Gene ID)	5013
蛋白 ID (SwissProt ID)	P32242.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000
蛋白分子量 (Molecular Weight)	Calculated MW: 37 kDa; Observed MW: 37 kDa

研究背景 (Background)

Probably plays a role in the development of the brain and the sense organs. Can bind to the BCD target sequence (BTS): 5'-

产品名称: OTX1 Rabbit Polyclonal Antibody
产品货号: APRab00018




TCTAATCCC-3'.

研究领域 (Research Area)

Neuroscience

图片 (Image Data)



Western blot analysis of Otx1 in mouse brain lysates using OTX1 antibody.

注意事项 (Note)

For research use only .