

产品名称: Hsc70 Rabbit Polyclonal Antibody
产品货号: APRab00037



产品概述 (Summary)

产品名称 (Production Name)	Hsc70 Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ICC/IF,FC,IP
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
纯化方式 (Purification)	Affinity Chromatography

免疫原信息 (Immunogen)

基因名 (Gene Name)	HSPA8
别名 (Alternative Names)	Heat shock cognate 71 kDa protein (Heat shock 70 kDa protein 8)
基因 ID (Gene ID)	3312
蛋白 ID (SwissProt ID)	P11142.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000,IHC 1:50-1:100,ICC/IF 1:50-1:200,FC 1:50-1:100,IP 1:20-1:50
蛋白分子量 (Molecular Weight)	Calculated MW: 71 kDa; Observed MW: 71 kDa

研究背景 (Background)

Acts as a repressor of transcriptional activation. Inhibits the transcriptional coactivator activity of CITED1 on Smad-

产品名称: Hsc70 Rabbit Polyclonal Antibody
产品货号: APRab00037




mediated transcription. Chaperone. Component of the PRP19-CDC5L complex that forms an integral part of the spliceosome and is required for activating pre-mRNA splicing.


研究领域 (Research Area)

Signal Transduction


图片 (Image Data)



Western blot analysis of Calreticulin in (1) HeLa lysates; (2) C6 lysates using Hsc70 antibody.



Immunohistochemistry analysis of paraffin-embedded Human breast cancer using Hsc70 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.



Immunofluorescence analysis of Hsc70 in MCF-7 using Hsc70 antibody.

产品名称: Hsc70 Rabbit Polyclonal Antibody
产品货号: APRab00037



注意事项 (Note)

For research use only .