

产品名称: GFAP Rabbit Polyclonal Antibody
产品货号: APRab00038



产品概述 (Summary)

产品名称 (Production Name)	GFAP Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ICC/IF,IP
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
纯化方式 (Purification)	Affinity Chromatography

免疫原信息 (Immunogen)

基因名 (Gene Name)	GFAP
别名 (Alternative Names)	GFAP; FLJ45472; cb345; ALXDRD
基因 ID (Gene ID)	2670
蛋白 ID (SwissProt ID)	P14136.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000,IHC 1:50-1:100,ICC/IF 1:50-1:200,IP 1:20-1:50
蛋白分子量 (Molecular Weight)	Calculated MW: 50 kDa; Observed MW: 50 kDa

研究背景 (Background)

GFAP is commonly used as a marker for intracranial and intraspinal tumors arising from astrocytes. In addition, GFAP

产品名称: GFAP Rabbit Polyclonal Antibody
产品货号: APRab00038




intermediate filaments are also present in nonmyelin-forming Schwann cells in the peripheral nervous system


研究领域 (Research Area)

Neuroscience


图片 (Image Data)



Western blot analysis of GFAP in rat brain lysates using GFAP antibody.



Immunohistochemistry analysis of paraffin-embedded Human glioma using GFAP antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.



Immunofluorescence analysis of GFAP in SH-SY5Y using GFAP antibody.

产品名称: GFAP Rabbit Polyclonal Antibody
产品货号: APRab00038



注意事项 (Note)

For research use only .