

产品名称: Rab9 Rabbit Polyclonal Antibody
产品货号: APRab00049



产品概述 (Summary)

产品名称 (Production Name)	Rab9 Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB, ICC/IF
种属反应性 (Reactivity)	Human, Mouse, Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Rabbit IgG in phosphate buffered saline, pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
纯化方式 (Purification)	Affinity Chromatography

免疫原信息 (Immunogen)

基因名 (Gene Name)	RAB9A
别名 (Alternative Names)	RAB9A; DmRab9; Ras related protein Rab 9A
基因 ID (Gene ID)	9367
蛋白 ID (SwissProt ID)	P51151.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000, ICC/IF 1:50-1:200
蛋白分子量 (Molecular Weight)	Calculated MW: 23 kDa; Observed MW: 23 kDa

研究背景 (Background)

Involved in the transport of proteins between the endosomes and the trans Golgi network. Interacts (preferentially in its

产品名称: Rab9 Rabbit Polyclonal Antibody
产品货号: APRab00049




GTP-bound form) with GCC2 (via its GRIP domain).

研究领域 (Research Area)

Tags & Cell Markers

图片 (Image Data)



Western blot analysis of Rab9 in HepG2 lysates using Rab9 antibody.

注意事项 (Note)

For research use only .