

产品名称: Hsp90 alpha/beta Rabbit Polyclonal Antibody
产品货号: AP Rab00057



产品概述 (Summary)

产品名称 (Production Name)	Hsp90 alpha/beta Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ICC/IF,FC,IP
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
纯化方式 (Purification)	Affinity Chromatography

免疫原信息 (Immunogen)

基因名 (Gene Name)	HSP90AA1/HSP90AB1
别名 (Alternative Names)	HSP90AA1; HSP90A; HSPC1; HSPCA; Heat shock protein HSP 90-alpha; Heat shock 86 kDa; HSP 86; HSP86; Renal carcinoma antigen NY-REN-38; HSP90AB1; HSP90B; HSPC2; HSPCB; Heat shock protein HSP 90-beta; HSP 90; Heat shock 84 kDa; HSP 84; HSP84
基因 ID (Gene ID)	3320/3326
蛋白 ID (SwissProt ID)	P07900/P08238.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000,IHC 1:50-1:100,ICC/IF 1:50-1:200,FC 1:50-1:100,IP 1:20-1:50
蛋白分子量 (Molecular Weight)	Calculated MW: 85 kDa; Observed MW: 90 kDa

产品名称: Hsp90 alpha/beta Rabbit Polyclonal Antibody
产品货号: APRab00057




研究背景 (Background)

Molecular chaperone that promotes the maturation, structural maintenance and proper regulation of specific target proteins involved for instance in cell cycle control and signal transduction. Undergoes a functional cycle that is linked to its ATPase activity. This cycle probably induces conformational changes in the client proteins, thereby causing their activation.


研究领域 (Research Area)

Signal Transduction

图片 (Image Data)




Western blot analysis of Hsp90 in (1) HeLa lysates; (2) Jurkat lysates; (3) RAW264.7 lysates; (4) NIH/3T3 lysates; (5) PC-12 lysates; (6) C6 lysates using Hsp90 alpha/beta antibody.



Immunohistochemistry analysis of paraffin-embedded Human liver carcinoma using Hsp90 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.

产品名称: Hsp90 alpha/beta Rabbit Polyclonal Antibody
产品货号: APRab00057



Immunofluorescence analysis of Hsp90 alpha/beta in Hela using Hsp90 antibody.

注意事项 (Note)

For research use only .