

产品名称: Lamin B Receptor Rabbit Polyclonal Antibody
产品货号: APRab00063



产品概述 (Summary)

产品名称 (Production Name)	Lamin B Receptor Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ICC/IF,IP
种属反应性 (Reactivity)	Human,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
纯化方式 (Purification)	Affinity Chromatography

免疫原信息 (Immunogen)

基因名 (Gene Name)	LBR
别名 (Alternative Names)	PHA; LMN2R; TDRD18; DHCR14B; LBR
基因 ID (Gene ID)	3930
蛋白 ID (SwissProt ID)	Q14739.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000,IHC 1:50-1:100,ICC/IF 1:50-1:200,IP 1:20-1:50
蛋白分子量 (Molecular Weight)	Calculated MW: 71 kDa; Observed MW: 71 kDa

研究背景 (Background)

The protein encoded by this gene belongs to the ERG4/ERG24 family. It localized in the nuclear envelope inner membrane

产品名称: Lamin B Receptor Rabbit Polyclonal Antibody
产品货号: APRab00063




and anchors the lamina and the heterochromatin to the membrane. It may mediate interaction between chromatin and lamin B. Mutations of this gene has been associated with autosomal recessive HEM/Greenberg skeletal dysplasia. Alternative splicing occurs at this locus and two transcript variants encoding the same protein have been identified.

研究领域 (Research Area)

Tags & Cell Markers

图片 (Image Data)



Western blot analysis of Lamin B Receptor in Jurkat lysates using Lamin B Receptor antibody.

注意事项 (Note)

For research use only .