

产品名称: DYNLL1 Rabbit Polyclonal Antibody
产品货号: APRab00064



产品概述 (Summary)

产品名称 (Production Name)	DYNLL1 Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ICC/IF,IP
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
纯化方式 (Purification)	Affinity Chromatography

免疫原信息 (Immunogen)

基因名 (Gene Name)	DYNLL1
别名 (Alternative Names)	LC8; PIN; DLC1; DLC8; LC8a; DNCL1; hdlc1; DNCLC1
基因 ID (Gene ID)	8655
蛋白 ID (SwissProt ID)	P63167.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000,IHC 1:50-1:100,ICC/IF 1:50-1:200,IP 1:20-1:50
蛋白分子量 (Molecular Weight)	Calculated MW: 10 kDa; Observed MW: 10 kDa

研究背景 (Background)

Acts as one of several non-catalytic accessory components of the cytoplasmic dynein 1 complex that are thought to be

产品名称: DYNLL1 Rabbit Polyclonal Antibody
产品货号: APRab00064




involved in linking dynein to cargos and to adapter proteins that regulate dynein function. Cytoplasmic dynein 1 acts as a motor for the intracellular retrograde motility of vesicles and organelles along microtubules. May play a role in changing or maintaining the spatial distribution of cytoskeletal structures.

研究领域 (Research Area)

Cell Biology

图片 (Image Data)



Western blot analysis of DYNLL1 in HeLa lysates using DYNLL1 antibody.

注意事项 (Note)

For research use only .