

产品名称: Cytochrome C Rabbit Polyclonal Antibody
产品货号: APRab00073



产品概述 (Summary)

产品名称 (Production Name)	Cytochrome C Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,IP
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
纯化方式 (Purification)	Affinity Chromatography

免疫原信息 (Immunogen)

基因名 (Gene Name)	CYCS
别名 (Alternative Names)	CYCS; CYC; Cytochrome c
基因 ID (Gene ID)	54205
蛋白 ID (SwissProt ID)	P99999.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000,IHC 1:50-1:100,IP 1:20
蛋白分子量 (Molecular Weight)	Calculated MW: 12 kDa; Observed MW: 15 kDa

研究背景 (Background)

CYCS Electron carrier protein. The oxidized form of the cytochrome c heme group can accept an electron from the heme

产品名称: Cytochrome C Rabbit Polyclonal Antibody
产品货号: APRab00073




group of the cytochrome c1 subunit of cytochrome reductase. Cytochrome c then transfers this electron to the cytochrome oxidase complex, the final protein carrier in the mitochondrial electron-transport chain.


研究领域 (Research Area)

Cardiovascular


图片 (Image Data)



Western blot analysis of Cytochrome C in HeLa lysates using Cytochrome C antibody.



Immunohistochemistry analysis of paraffin-embedded Human colonic cancer using Cytochrome C antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.



Immunofluorescence analysis of Cytochrome C in HeLa cells using Cytochrome C antibody.

产品名称: Cytochrome C Rabbit Polyclonal Antibody
产品货号: APRab00073



注意事项 (Note)

For research use only .