

产品名称: Cytochrome C Rabbit Polyclonal Antibody
产品货号: APRab00095



产品概述 (Summary)

产品名称 (Production Name)	Cytochrome C Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ICC/IF,IP
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
纯化方式 (Purification)	Affinity Chromatography

免疫原信息 (Immunogen)

基因名 (Gene Name)	CYCS
别名 (Alternative Names)	CYCS; CYC; Cytochrome c
基因 ID (Gene ID)	54205
蛋白 ID (SwissProt ID)	P99999.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000,IHC 1:50-1:100,ICC/IF 1:50-1:200,IP 1:20-1:50
蛋白分子量 (Molecular Weight)	Calculated MW: 12 kDa; Observed MW: 12 kDa

研究背景 (Background)

CYCS Electron carrier protein. The oxidized form of the cytochrome c heme group can accept an electron from the heme

产品名称: Cytochrome C Rabbit Polyclonal Antibody
产品货号: APRab00095




group of the cytochrome c1 subunit of cytochrome reductase. Cytochrome c then transfers this electron to the cytochrome oxidase complex, the final protein carrier in the mitochondrial electron-transport chain.


研究领域 (Research Area)

Cardiovascular


图片 (Image Data)



Western blot analysis of Cytochrome C in Human kidney lysates using Cytochrome C antibody.



Immunohistochemistry analysis of paraffin-embedded Human liver using Cytochrome C antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.



Immunofluorescence analysis of Cytochrome C in HeLa using Cytochrome C antibody.

产品名称: Cytochrome C Rabbit Polyclonal Antibody
产品货号: APRab00095



注意事项 (Note)

For research use only .