

产品名称: Rel B Rabbit Polyclonal Antibody
产品货号: APRab00099



产品概述 (Summary)

产品名称 (Production Name)	Rel B Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ICC/IF,IP
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
纯化方式 (Purification)	Affinity Chromatography

免疫原信息 (Immunogen)

基因名 (Gene Name)	RELB
别名 (Alternative Names)	RELB; Transcription factor RelB; I-Rel
基因 ID (Gene ID)	5971
蛋白 ID (SwissProt ID)	Q01201.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000,IHC 1:50-1:100,ICC/IF 1:50-1:200,IP 1:20-1:50
蛋白分子量 (Molecular Weight)	Calculated MW: 62 kDa; Observed MW: 70 kDa

研究背景 (Background)

Transcription factors of the nuclear factor κ B (NF- κ B)/Rel family play a pivotal role in inflammatory and immune responses.

产品名称: Rel B Rabbit Polyclonal Antibody
产品货号: APRab00099




Involved in many biological processes such as inflammation, immunity, differentiation, cell growth, tumorigenesis and apoptosis.


研究领域 (Research Area)

Cell Biology

图片 (Image Data)



Western blot analysis of RelB in Raji lysates using Rel B antibody.



Immunohistochemistry analysis of paraffin-embedded Human tonsil using RelB antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.

注意事项 (Note)

For research use only .