

产品名称: INPP4B Rabbit Polyclonal Antibody
产品货号: APRab00101



产品概述 (Summary)

产品名称 (Production Name)	INPP4B Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,FC
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
纯化方式 (Purification)	Affinity Chromatography

免疫原信息 (Immunogen)

基因名 (Gene Name)	INPP4B 4-bisphosphate 4-phosphatase; Inositol polyphosphate 4 phosphatase type II
别名 (Alternative Names)	105kDa; Inositol polyphosphate 4-phosphatase type II; Inpp4b; MGC132014; Type II inositol 3; 4 bisphosphate 4 phosphatase; Type II inositol-3
基因 ID (Gene ID)	8821
蛋白 ID (SwissProt ID)	O15327.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000,IHC 1:50-1:100,FC 1:50-1:100
蛋白分子量 (Molecular Weight)	Calculated MW: 105 kDa; Observed MW: 105 kDa

产品名称: INPP4B Rabbit Polyclonal Antibody
产品货号: APRab00101




研究背景 (Background)

Catalyzes the hydrolysis of the 4-position phosphate of phosphatidylinositol 3,4-bisphosphate, inositol 1,3,4-trisphosphate and inositol 1,4-bisphosphate.

研究领域 (Research Area)

Signal Transduction

图片 (Image Data)



Western blot analysis of INPP4B in MCF-7 lysates using INPP4B antibody.

注意事项 (Note)

For research use only .