

产品名称: IRF3 Rabbit Polyclonal Antibody
产品货号: APRab00102



产品概述 (Summary)

产品名称 (Production Name)	IRF3 Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ICC/IF,FC
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
纯化方式 (Purification)	Affinity Chromatography

免疫原信息 (Immunogen)

基因名 (Gene Name)	IRF3
别名 (Alternative Names)	IRF3; Interferon regulatory factor 3; IRF-3
基因 ID (Gene ID)	3661
蛋白 ID (SwissProt ID)	Q14653.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000,IHC 1:50-1:100,ICC/IF 1:50-1:200,FC 1:50-1:100
蛋白分子量 (Molecular Weight)	Calculated MW: 47 kDa; Observed MW: 57 kDa

研究背景 (Background)

IRFs comprise a family of transcription factors that function within the Jak/Stat pathway to regulate interferon (IFN) and

产品名称: IRF3 Rabbit Polyclonal Antibody
产品货号: APRab00102




IFN-inducible gene expression in response to viral infection. IRF-3 can inhibit cell growth and plays a critical role in controlling the expression of genes in the innate immune response.


研究领域 (Research Area)

Signal Transduction


图片 (Image Data)



Western blot analysis of IRF3 in (1) HeLa lysates; (2) 3T3 lysates using IRF3 antibody.



Immunohistochemistry analysis of paraffin-embedded Human cervix carcinoma using IRF3 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.



Immunofluorescence analysis of IRF3 in Jurkat using IRF3 antibody.

产品名称: IRF3 Rabbit Polyclonal Antibody
产品货号: APRab00102



注意事项 (Note)

For research use only .