

产品名称: IL-1 Receptor Antagonist Protein Rabbit Polyclonal Antibody
产品货号: APRab00104



产品概述 (Summary)

产品名称 (Production Name)	IL-1 Receptor Antagonist Protein Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ICC/IF,FC,IP
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
纯化方式 (Purification)	Affinity Chromatography

免疫原信息 (Immunogen)

基因名 (Gene Name)	IL1RN
别名 (Alternative Names)	IL-1 Receptor Antagonist Protein; DIRA; IRAP; IL1F3; IL1RA; MVCD4; IL-1RN; IL-1ra; IL-1ra3; ICIL-1RA
基因 ID (Gene ID)	3557
蛋白 ID (SwissProt ID)	P18510.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000,IHC 1:50-1:100,ICC/IF 1:50-1:200,FC 1:50-1:100,IP 1:20-1:50
蛋白分子量 (Molecular Weight)	Calculated MW: 20 kDa; Observed MW: 20 kDa

**产品名称: IL-1 Receptor Antagonist Protein Rabbit
Polyclonal Antibody**
产品货号: APRab00104




研究背景 (Background)

Interleukin-1 (IL-1) is a key pro-inflammatory cytokine that has diverse actions in the brain as a regulator of host defense responses and a mediator of inflammation. Inhibits the activity of interleukin-1 by binding to receptor IL1R1 and preventing its association with the coreceptor IL1RAP for signaling. Has no interleukin-1 like activity.

研究领域 (Research Area)

Immunology

图片 (Image Data)



Western blot analysis of IL-1RA in (1) HeLa lysates; (2) 3T3 lysates using IL-1 Receptor Antagonist Protein antibody.

注意事项 (Note)

For research use only .