

产品名称: ARF1 Rabbit Polyclonal Antibody
产品货号: APRab00105



产品概述 (Summary)

产品名称 (Production Name)	ARF1 Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,FC
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
纯化方式 (Purification)	Affinity Chromatography

免疫原信息 (Immunogen)

基因名 (Gene Name)	ARF1
别名 (Alternative Names)	ADP-ribosylation factor 1
基因 ID (Gene ID)	375
蛋白 ID (SwissProt ID)	P84077.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000,IHC 1:50-1:100,FC 1:50-1:100
蛋白分子量 (Molecular Weight)	Calculated MW: 21 kDa; Observed MW: 21 kDa

研究背景 (Background)

Stimulate the ADP-ribosyltransferase activity of cholera toxin and play a role in vesicular trafficking as activators of

产品名称: ARF1 Rabbit Polyclonal Antibody
产品货号: APRab00105




phospholipase D. GTP-binding protein that functions as an allosteric activator of the cholera toxin catalytic subunit, an ADP-ribosyltransferase. Involved in protein trafficking among different compartments.

研究领域 (Research Area)

Signal Transduction

图片 (Image Data)



Western blot analysis of ARF1 in HEK293 lysates using ARF1 antibody.

注意事项 (Note)

For research use only .